

PROSTHETICS



DENTAL CATALOG

WHAT DO YOU CALL A CONVICTION
THAT EVERYTHING CAN ALWAYS BE
A LITTLE BETTER, EASIER AND MORE
RELIABLE?

WE CALL IT
Upgrade Dentistry

With expertise and enthusiasm, COLTENE turned this claim into reality. We know exactly our end customers needs and want and work tirelessly to provide them with innovative and smart solutions.

This is not an empty promise, it is our mission. Upgrading Dentistry makes dental treatments more effective, safe and economical for dental professionals and patients.

Better Quality. Better Reliability. Better Practise.



www.coltene.com

1 RETRACTION

ROEKO Comprecap / Comprecap anatomic	10
ROEKO Stay-put	10
ROEKO Retracto	11
Retraction Accessories	11

2 BITE REGISTRATION

JET BLUE BITE	14
JET BITE	14

3 IMPRESSION A-SILICONE

AFFINIS	18
AFFINIS Light Body / Fast	19
AFFINIS Regular Body / Fast	19
AFFINIS PRECIOUS	19
AFFINIS Heavy Body / Fast	20
AFFINIS MonoBody	20
AFFINIS BLACK EDITION	21
AFFINIS System 360 Putty	21
AFFINIS Putties	21
PRESIDENT The Original	22
PRESIDENT The Original Xtra Light Body	23
PRESIDENT The Original Light Body	23
PRESIDENT The Original Regular Body	23
PRESIDENT The Original System 360 MonoBody	23
PRESIDENT The Original Heavy Body	24
PRESIDENT The Original Putties	24

4 IMPRESSION C-SILICONE

Speedex	28
Speedex Light Body	29
Speedex Medium	29

Speedex Putty	29
Speedex Putty Soft	29
Speedex Activator	30
Pressure Spot Indicator	30

5 PROSTHETICS ACCESSORIES

COLTENE Adhesive	34
Impression Trays	35
Impression Trays Sets	35
Impression Trays Quadrant & Anterior	36
Mixing Devices	36
Mixing Tips	36
Other Accessories	37

6 LABORATORY

GI-MASK	40
GI-MASK Automix New Formula	40
GI-Mask 3D	41
GI-MASK Accessories	41
Lab Putty	42
Lab Putty / Lab Putty Hard	43
ACCU-TRAC	44
ACCU-TRAC	45
PINDEX	46
PINDEX Laser System	46
Drills	47
Pins & Accessories	47
Saw Blades	47
HYGENIC Waxes	48
HYGENIC bite-wafers	48
HYGENIC u-shaped occlusal rim wax	49
HYGENIC yellow bite wax	49
HYGENIC boxing wax	49
HYGENIC medium-soft wax	49
HYGENIC orthodontic tray wax	50
HYGENIC wax ropes	50
HYGENIC wax sticks	50

TABLE OF CONTENTS

HYGENIC extra tough wax	50
Laboratory Accessories	51
Laboratory Accessories	52

SPOTLIGHT



ROEKO Comprecap

10

ROEKO Comprecap compression caps stop bleeding naturally by compression and control the moisture around the tooth stump.

ROEKO Comprecap / Comprecap anatomic



530000



530000

ROEKO Comprecap compression caps stop bleeding naturally by compression and control the moisture. The use is simple and supports impression preparation regardless of the retraction technique. ROEKO Comprecap anatomic are formed compression caps with semicircled shaped spaces on two opposite sides and can therefore be easily placed on adjacent teeth.

REF	Item
530000	Set, 150 pcs
530008	Size 1, 120 pcs
530010	Size 3, 120 pcs
530014	Size 5, 60 pcs
531008	Size 1, 600 pcs
531010	Size 3, 600 pcs
531014	Size 5, 400 pcs
531201	Anatomic Size 1, 120 pcs
531203	Anatomic Size 3, 120 pcs
531205	Anatomic Size 5, 60 pcs
531299	Anatomic Set, 120 pcs

ROEKO Stay-put



521000

Stay-put combines the advantages of a braided retraction cord with the adaptability of a fine copper filament. The pliable core is so effective that the cord is not only easy to place in the sulcus but it stays there.

REF	Item
521000	Non-Impregnated, X-Fine, 1 pcs
521001	Non-Impregnated, Fine, 1 pcs
521002	Non-Impregnated, Medium, 1 pcs
521003	Non-Impregnated, Thick, 1 pcs

RETRACTION

ROEKO Retracto

Retracto cords are available as non impregnated. The mechanical use of the retraction cord causes compression. The compression contracts the treated tissue and causes haemostasis.

Item	REF				
	Braided		Twisted		
ROEKO Retracto	Non Impregnated	x-fine	1 pcs	511000	523000
ROEKO Retracto	Non Impregnated	fine	1 pcs	511001	523001
ROEKO Retracto	Non Impregnated	medium	1 pcs	511002	523002
ROEKO Retracto	Non Impregnated	thick	1 pcs	511003	-



511000



523000

Retraction Accessories

Double ended, stainless steel instruments facilitate placing of retraction cord around the tooth.

- Serrated or smooth edges
- Structured handle
- Ergonomically designed to work around the tooth
- The balanced notches hold the cord in place
- Sterilizable



520000

REF	Item
520000	Retraction Cord Packer, Smooth, 1 pcs
520001	Retraction Cord Packer, Serrated, 1 pcs



2

Bite Registration

SPOTLIGHT



JET BLUE BITE

14

Jet Blue Bite is a material for bite registration. It offers a fast intra oral setting-time and sets rock hard for greater accuracy but is still flexible enough for easy trimming or cutting.

JET BITE

14

Jet Bite is a material for bite registration. The direct placement and quick intra oral setting time as well as the ideal stability of the material make it easy to record the patients occlusion accurately.

JET BLUE BITE



Jet Blue Bite bite registration material flows with a smooth whipped-cream consistency exactly where you want it in a predictable quantity and ideal ribbon pattern. It offers a fast oral setting-time and sets rock hard for greater accuracy but is still flexible enough for easy trimming or cutting. The superior handling characteristics are the key for a perfect occlusal fit on the first try.

- Stays exactly where you want it
- No resistance to closure, no rebound!
- Rock hard for accuracy, yet flexible for easy trimming



REF	Item	
6491	Fast Refill	
	1 × 50 ml	Cartridge
	1 × 6 pcs	Universal Mixing Tip
6493	Superfast Refill	
	1 × 50 ml	Cartridge
	1 × 6 pcs	Universal Mixing Tip
6495	Fast Refill	
	4 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip
6497	Superfast Refill	
	4 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip

JET BITE



Jet Bite gives a great occlusal fit on the first try. The patient bites down naturally every time because of Jet Bite’s fluffy, mousse-like consistency. There is no need for time-consuming adjustments. Jet Bite flows smoothly - sets rock hard.

- Smooth mousse-like consistency - No resistance to closure, no rebound
- High dimensional stability - no distortion
- Not brittle - easy to grind and cut

REF	Item	
6400	Refill	
	2 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip
6410	Refill, 20 × 50 ml	



NOTES



3

Impression A-Silicone

SPOTLIGHT



AFFINIS

18

AFFINIS is COLTENEs premium A-Silicone for dental impressions. Apart from the ideal working and setting time, AFFINIS convinces dentists with a high dimensional stability and accuracy. The extraordinary surfactant system results in outstanding flow properties for greater detail reproduction.

PRESIDENT The Original

22

The launch of PRESIDENT, the first A-silicone on the market, is considered a milestone for dental impression materials. Due to the attractive portfolio width, PRESIDENT The Original offers a perfect choice of consistencies with excellent flow properties in moist conditions.

AFFINIS

Perfect Impressions - The crowning in level of detail

The perfect AFFINIS impression is the key to accurately fitting dental restorations. It forms the basis for successful teamwork between dentist and dental technician.

Still one of the fastest A-silicones

AFFINIS fast, with a maximum total setting time of 2 min 10 sec guarantees optimum efficiency and saves time.

Precise moulding of every detail even under the most difficult conditions

The highly active surfactant system is characterized by excellent surface affinity in a moist environment.

Confidence during impression taking

Self-contouring consistency prevents air bubbles and distortions.

Easy monitoring for success

AFFINIS PRECIOUS wash materials with impressive gold and silver pigmentation for excellent readability and detail reproduction.

Autoclavable impressions

All AFFINIS and AFFINIS PRECIOUS impressions can be autoclaved at 134°C without losing its unique properties.



AFFINIS Light Body / Fast

6501

AFFINIS Light Body is an extremely light-flowing correction material that can easily capture the smallest details.

REF	Item
6501	Refill Light Body 2 × 50 ml Cartridge 12 pcs Universal Mixing Tip
6580	Refill Light Body, 20 × 50 ml
6601	Refill Fast Light Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6680	Refill Fast Light Body, 20 × 50 ml

AFFINIS Regular Body / Fast

6511

AFFINIS Regular Body is a medium-flowing correction material with good stability.

REF	Item
6511	Refill Regular Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6582	Refill Regular Body, 20 × 50 ml
6611	Refill Fast Regular Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6682	Refill Fast Regular Body, 20 × 50 ml

AFFINIS PRECIOUS

6775

The AFFINIS PRECIOUS Light and Regular Body convinces with its coordinated flowability and features a very low contact angle. The unique silver and gold pigmentation significantly reduces light diffusion and improves the visual impression of details.

REF	Item
6775	Light Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6776	Regular Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip

AFFINIS Heavy Body / Fast



6590

AFFINIS Heavy Body is a viscous tray material, available in 50ml and 75ml cartridges as well as in 5:1 cartridges for use in the Automix system. With its surface affinity and optimal stability during rapid pressure build-up, it is the perfect tray material.

REF	Item
6590	Refill Heavy Body 2 × 50 ml Cartridge 1 × 6 pcs Mixing Tips
6595	Refill Heavy Body, 20 × 50 ml
6520	Refill Heavy Body 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
6585	Refill Heavy Body, 20 × 75 ml
6485	System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
6487	System 360 Refill, 2 × 380 ml
6620	Refill Fast Heavy Body 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
6650	Refill Fast Heavy Body 2 × 50 ml Cartridge 1 × 6 pcs Mixing Tips
6652	Refill Fast Heavy Body, 20 × 50 ml
6685	Refill Fast Heavy Body, 20 × 75 ml

AFFINIS MonoBody



6480

AFFINIS MonoBody is outstanding for its high stability and thixotropy. Due to its adapted flow properties, it is ideal for implant and/or single-phase impressions.

REF	Item
6480	System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
6482	System 360 Refill, 2 × 380 ml
6547	Refill 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
6587	Refill, 20 × 75 ml

AFFINIS BLACK EDITION



60019554

The black tray material AFFINIS Heavy Body BLACK EDITION contrasts the gold or silver colored AFFINIS PRECIOUS and ensures excellent viewability of the impression.

REF	Item	
60019554	Refill	
	2 × 75 ml	Cartridge
60019774	1 × 8 pcs	Mixing Tips
	System 360 Starter Kit	
60019776	1 × 380 ml	Cartridge
	1 × 10 pcs	Dynamic Mixing Tips
60019776	System 360 Refill, 2 × 380 ml	

AFFINIS System 360 Putty



6471

AFFINIS System 360 putty has a strong “firm body” which is immediately noticeable once dispensed from the mixing unit with no signs of stickiness. The dentist is comforted with a reliable, automatically-mixed putty material for successful precision impressions.

REF	Item	
6471	Putty System 360 Starter Kit	
	1 × 380 ml	Cartridge
6472	1 × 10 pcs	Dynamic Mixing Tips
	System 360 Refill, 2 × 380 ml	

AFFINIS Putties



6530

The AFFINIS putties are highly viscous, kneadable tray materials. The remarkable recoverability of Putty Soft and Putty Super Soft leads to highly precise impressions even in areas with strong undercuts. The compatible consistencies of the base and catalyst can be mixed quickly and easily into a homogeneous material.

REF	Item	
6530	Putty Soft	
	1 × 300 ml	Base
6531	1 × 300 ml	Catalyst
	2 ×	Dosage spoon
6531	Fast Putty Soft	
	1 × 300 ml	Base
6535	1 × 300 ml	Catalyst
	2 ×	Dosage spoon
6535	Putty Super Soft	
	1 × 300 ml	Base
6535	1 × 300 ml	Catalyst
	2 ×	Dosage spoon

PRESIDENT The Original

A small step for dentists. A great step for dentistry.

With the launch of PRESIDENT in 1975 - the first A-silicone on the market - a milestone in the development of high-precision dental impression materials was created. Through continuous further development, the PRESIDENT The Original brand continues to impress with its remarkable quality today.

The highly pigmented colors of the correction materials provide better contrast and show every detail, no matter how small, even before the impression is poured. PRESIDENT The Original offers a perfect range of consistencies with excellent flow properties in moist environments.

- High tensile strength and elasticity for more safety when taking impressions
- Excellent flow properties and thixotropy for impression taking of even the smallest details
- High-contrast colors for optimal viewability



PRESIDENT The Original Xtra Light Body



60019937

PRESIDENT the Original Xtra Light Body is a light-flowing material that can easily capture even the smallest marginal details.

REF	Item	
60019937	Refill	
	2 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip

PRESIDENT The Original Light Body



60019938

PRESIDENT the Original Light body is characterized by its excellent flow properties and viscosity. As a real all-rounder it is suitable for any impression technique.

REF	Item	
60019938	Refill	
	2 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip
60019949	Tube Refill	
	1 × 53 ml	Tube Base
	1 × 53 ml	Tube Catalyst

PRESIDENT The Original Regular Body



60019939

PRESIDENT the Original Regular Body is a medium-flowing material which is very suitable for the double mixing technique due to its stability.

REF	Item	
60019939	Refill	
	2 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip
60019950	Tube Refill	
	1 × 53 ml	Tube Base
	1 × 53 ml	Tube Catalyst

PRESIDENT The Original System 360 MonoBody



60019945

PRESIDENT the Original MonoBody is ideally suited for implant and/or single-phase impressions due to its adapted flow properties. It is available in 5:1 cartridges for use in the Automix System.

REF	Item	
60019945	System 360 Starter Kit	
	1 × 380 ml	Cartridge
	1 × 10 pcs	Dynamic Mixing Tips
60019946	System 360 Refill, 2 × 380 ml	

PRESIDENT The Original Heavy Body



60019940

PRESIDENT the Original Heavy Body is a viscous tray material available in 75ml cartridges for use with a dispenser and 5:1 cartridges for use in the Automix System.

REF	Item	
60019940	Refill	
	2 × 50 ml	Cartridge
	1 × 6 pcs	Mixing Tips
60019941	Refill	
	2 × 75 ml	Cartridge
	1 × 8 pcs	Mixing Tips
60019947	System 360 Starter Kit	
	1 × 380 ml	Cartridge
	1 × 10 pcs	Dynamic Mixing Tips
60019948	System 360 Refill, 2 × 380 ml	

PRESIDENT The Original Putties



60019942



60019943

PRESIDENT the Original Putties are kneadable, highly viscous tray materials. The remarkable recoverability of Putty Soft and Putty Super Soft results in highly precise impressions even in areas with strong undercuts.

REF	Item	
60019942	Putty	
	1 × 300 ml	Base
	1 × 300 ml	Catalyst
	2 ×	Dosage spoon
60019943	Putty Soft	
	1 × 300 ml	Base
	1 × 300 ml	Catalyst
	2 ×	Dosage spoon
60022166	Fast Putty Soft	
	1 × 300 ml	Base
	1 × 300 ml	Catalyst
	2 ×	Dosage spoon
60019944	Putty Super Soft	
	1 × 300 ml	Base
	1 × 300 ml	Catalyst
	2 ×	Dosage spoon



NOTES



4

Impression C-Silicone

SPOTLIGHT



Speedex

28

Speedex is a premium C-silicone with an A-silicone quality. Due to the low shrinkage, the impressions remain dimensionally stable and can be poured without problems for up to 7 days.

Speedex

Bites of passion

Speedex has elastic recoverability and reproduction, which is otherwise only achieved by A-Silicones. As a result of the low shrinkage, impressions remain dimensionally stable and can be cast for up to 7 days without problems. Furthermore this premium C-Silicone shows an outstanding wettability and excellent flow characteristics. Saliva and blood are displaced by Speedex medium and have no negative side effects.

The Speedex system gives the advantage to work with only one activator for 4 consistencies.

- 7 days dimension stability
- Controllable working time
- Surface active
- Moisture control
- Optimized precision



Speedex Light Body



4980

Speedex Light Body is outstanding for its good flow properties and accurately molds every detail. The thin-flowing consistency is mainly used as a correction compound for the correction impression technique and as an injection compound in the two-phase impression technique.

REF	Item
4980	Light Body, 1 × 140 ml

Speedex Medium



4986

Speedex Medium is a medium-flow material that is ideal for single-phase impressions and as a relining impression material due to its flow properties.

REF	Item
4986	Medium, 1 × 140 ml

Speedex Putty



4970

Speedex Putty is a kneadable tray material which is suitable for all kinds of impression techniques.

REF	Item
4970	Putty, 1 × 910 ml

Speedex Putty Soft



60023572

Speedex Putty Soft is smooth and easy to knead to a homogeneous consistency.

REF	Item
60023572	Putty soft, 1 × 910 ml

Speedex Activator

Universal Activator for Speedex Putty, Speedex Putty Soft, Speedex Light Body und Speedex Medium.



4990

REF	Item
4990	Universal Activator, 1 × 60 ml

Pressure Spot Indicator

The polysiloxane-based material for the location of high spots in dentures.



4440

- Recognition of pressure points on new or existing
- partial and total dentures
- Determining the need for relining
- Testing of the exact fitting of inlays, crowns, bridges and cast brackets

REF	Item
4440	PSI (Base & Catalyst) 1 × 55 ml Tube Base 1 × 20 ml Tube Catalyst
4445	PSI Activator, 1 × 20 ml

5

Prosthetics Accessories

SPOTLIGHT



Impression Trays

35

PRESIDENT Trays are single use plastic impression trays. A patented system of intersecting ribs and rails provides unprecedented rigidity in a plastic tray. The PRESIDENT Trays AC are also autoclavable.

Mixing Tips

36

Universal Mixing Tip - a universal mixer developed especially for all light and regular body materials as well as bite registration and gingival mask materials from COLTENE. The reduced extrusion force allows for a comfortable and accurate dispensing of material.

COLTENE Adhesive



4405

A thin layer of COLTENE Adhesive ensures adhesion of tray material to the tray sides.

COLTENE Adhesive

COLTENE Adhesive is a universal adhesive for silicone-based dental impression materials and serves as an adhesive bonding agent between the impression material and the impression tray. Separation of the impression material during removal from the mouth can therefore be prevented.

COLTENE Adhesive AC

COLTENE Adhesive is an autoclavable adhesive for all AFFINIS impression materials and serves as an adhesive bonding agent between this material and the PRESIDENT disposable impression trays AC. Separation of the impression material during removal from the mouth can therefore be prevented. The special spray nozzle provides application with a uniform, thin layer without becoming clogged. COLTENE Adhesive AC was specially developed for the autoclavable system (AFFINIS impression materials, PRESIDENT disposable impression trays AC).

REF	Item
4405	Adhesive Autoclavable, 1 × 10 ml
4410	Adhesive, 1 × 10 ml

Impression Trays

A patented system of intersecting ribs and rails provides unprecedented rigidity in a plastic tray. The anatomic form relieves the fitting and ensures a perfect impression. Autoclavable and non autoclavable.

- Anatomically correct: A statistical database accumulated over 20 years has enabled COLTENE to define the ideal form for each tray system
- Built-in retention: The PRESIDENT Tray System's design aggressively retains all types of impression material.
- Transparent: Easy and precise positioning of the tray.
- Modifiable form: The tray can individually be adjusted with heat.
- Single use: No risk of contamination, no time-consuming cleaning.

Item			REF					
			Upper Jaw	Upper Jaw	Upper Jaw	Lower Jaw	Lower Jaw	Lower Jaw
			S	M	L	S	M	L
Impression Trays	12 pcs	Autoclavable	C6954	C6952	C6950	C6955	C6953	C6951
		Not autoclavable	C6905	C6903	C6901	C6906	C6904	C6902



C6950



C6901

Impression Trays Sets



C6955

REF	Item
C6910	Assortment Impression Trays 1 × 12 pcs Full Upper Jaw Large 1 × 12 pcs Full Lower Jaw Large 1 × 12 pcs Full Upper Jaw Medium 1 × 12 pcs Full Lower Jaw Medium 1 × 12 pcs Full Upper Jaw Small 1 × 12 pcs Full Lower Jaw Small 1 × Quadrant UJ le. / LJ ri. 1 × Quadrant UJ ri. / LJ le. 1 × Anterior UJ / LJ
C6959	Assortment PRESIDENT Trays, autoclavable 1 × 12 pcs Full Upper Jaw Large, autoclavable 1 × 12 pcs Full Lower Jaw Large, autoclavable 1 × 12 pcs Full Upper Jaw Medium, autoclavable 1 × 12 pcs Full Lower Jaw Medium, autoclavable 1 × 12 pcs Full Upper Jaw Small, autoclavable 1 × 12 pcs Full Lower Jaw Small, autoclavable 1 × Quadrant UJ le. / LJ ri., autoclavable 1 × Quadrant UJ ri. / LJ le., autoclavable 1 × Anterior UJ / LJ, autoclavable

Impression Trays Quadrant & Anterior



C6956

REF	Item
C6907	Quadrant UJ le. / LJ ri., 12 pcs
C6908	Quadrant UJ ri. / LJ le., 12 pcs
C6909	Anterior UJ / LJ, 12 pcs
C6956	Quadrant UJ le. / LJ ri., autoclavable, 12 pcs
C6957	Quadrant UJ ri. / LJ le., autoclavable, 12 pcs
C6958	Anterior UJ / LJ, autoclavable, 12 pcs

Mixing Devices



6985

Dispensers and other mixing devices for cartridges

REF	Item
6985	Plunger 50/75 ml, 1 pcs
60011147	Sympress II Europe, 1 pcs
60011150	Sympress II UK, 1 pcs

REF	Accessories
6900	Dispenser 50/75 ml, 1 pcs

Mixing Tips



6552

COLTENE Mixing Tips and Oral Tips for different wash and tray materials which ensure a homogenous ratio between base and catalyst.

REF	Item
6162	Dynamic Mixing Tip Yellow, 40 pcs
6163	System 360 Fixation Ring yellow, 1 pcs
6210	Oral Tip white, 100 pcs
6223	Spreader Tip, 100 pcs
6552	Mixing Tip turquoise, 40 pcs
6555	Oral Tip yellow, 100 pcs
60019663	Universal Mixing Tip Refill, 40 pcs
60019694	Universal Mixing Tip Refill, 100 pcs

Other Accessories



65002689

REF	Item
4420	Cutting Knife, 1 pcs
4412	Mixing Pad, 1 pcs
4414	Mixing Pad, 10 pcs
65002689	Coltopress, 1 pcs



6

Laboratory

SPOTLIGHT



GI-MASK

41

GI-MASK is a silicone-based gingival mask with which the exact course of the gingiva is recorded in the preparation area and remains secured until the denture is completed. The aesthetic coloring of the materials helps you to better recreate the clinical situation.

Lab Putty

42

Lab-Putty is the outstanding polysiloxane laboratory modelling material with a 62 - 67 Shore A hardness. Ideal for model repairs, isolation, denture heat protection, keys, block outs and matrixes for the formation of general occlusion and functional guides.

GI-MASK



8030



8036

GI-MASK is a C-Silicone based gingival mask for precise and fast reproduction of gingival contours. The cured material is dimensionally stable for several days. Another advantage of GI-MASK is the exclusive application system, which stands for safe dosing, clean mixing by hand and easy application.

REF	Item	
8030	Starter Kit	
	1 × 150 ml	Base
	1 × 18 ml	Activator
	1 × 50 ml	Universal Separator
8036	Activator, 1 × 18 ml	
8037	Refill	
	1 × 150 ml	Base
	1 × 18 ml	Activator

GI-MASK Automix New Formula



8060

GI-MASK Automix New Formula is a gingival mask based on A-Silicone which is used for precise detection of the gingival line. Due to the practical cartridge format, it is time-saving and clean to work with.

REF	Item	
8060	Starter Kit	
	2 × 50 ml	Cartridge
	1 × 50 ml	Universal Separator
	1 ×	Trimmer for Silicones
	1 × 12 pcs	Universal Mixing Tip
	1 × 12 pcs	Oral Tips
8065	Refill	
	2 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip
	1 × 12 pcs	Oral Tips

GI-Mask 3D



60022251

With the scannable gingival mask material GI-MASK 3D, reproduced gingival contours in the preparation area can be digitized without scan spray. The use of a gingival mask is particularly important in the fabrication of implant-supported superstructures in order to obtain the most accurate-fitting dental prosthesis possible.

REF	Item	
60022250	STARTER KIT	
	2 × 50 ml	Cartridge
	1 × 50 ml	Universal Separator
	1 ×	Trimmer for Silicones
	1 × 12 pcs	Universal Mixing Tip
	1 × 12 pcs	Oral Tips
60022251	Refill	
	2 × 50 ml	Cartridge
	1 × 12 pcs	Universal Mixing Tip
	1 × 12 pcs	Oral Tips

GI-MASK Accessories



8000

Accessories for GI-MASK

Universal Separator, separating liquid for addition- and condensation type silicones in dental laboratories

REF	Item
8000	Universal Separator, 1 × 50 ml

Lab Putty

COLTENEs Lab-Putty was especially designed for the use in laboratory. It can be used for model repairs, isolation, denture heat protection, keys, block outs and matrixes for formation of general occlusion and functional patterns. Lab-Putty hard is characterized by its high heat-resistance (300°C) and its high end hardness (90 Shore A).



Lab Putty / Lab Putty Hard



8805



8840

Lab-Putty is the outstanding polysiloxane laboratory modelling product material with 63 - 70 Shore A for model repairs, isolation and denture heat protection, keys and block outs and matrixes for formation of general occlusion and functional patterns.

A smooth acrylic surface is guaranteed with the use of Lab Putty hard for relining, keys and repair models. Separate isolation is eliminated. Its high heat-resistance (300° C) and perfect dimensional stability makes it suitable for soldering models (90 Shore A).

REF	Item	
8805	Lab Putty	
	1 × 900 ml	Base
	1 × 40 ml	Activator
	1 ×	Dosage spoon
8815	Lab Putty	
	1 × 5000 ml	Base
	5 ×	Dosage spoon
8880	Lab Putty Hard	
	1 × 2600 ml	Base
	5 ×	Dosage spoon
8882	Lab Putty Hard	
	1 × 900 ml	Base
	1 × 40 ml	Activator
	1 ×	Dosage spoon
8840	Activator, 1 × 40 ml	

ACCU-TRAC

The ACCU-TRAC system makes the fabrication of precision model and dies easier. It creates a precise base for your prosthodontic work. The model and base are poured as one, eliminating the additional work when using pins and sleeves. Time and material can be saved without sacrificing precision.



ACCU-TRAC



ATK25



ATA20

The ACCU-TRAC System makes the fabrication of precision models and dies easier: The model and base are poured as one, eliminating the additional work when using pins and sleeves. You save time and material, without sacrificing precision.

- One universal tray size, accomodates all maxillary and mandibular sizes, including quadrants
- Low base height, model fits any articulator
- One pour technique saves time and materials
- Economical, all components reusable
- Numerical die location, immediate recognition of dies, easy die repositioning

REF	Item	
AT6	Introductory Kit	
	1 × 6 pcs	Precision Die Trays
	1 × 6 pcs	Pouring/Ejector Bases
	1 × 6 pcs	Articulation Spacers
	1 × 25 pcs	Articulation Keepers
AT72	Economy Pack, 72 pcs	
ATR12	Acrylic frame Refill	
	1 × 12 pcs	Precision Die Trays
	1 × 50 pcs	Articulation Keepers
ATR144	Economy Pack, 144 pcs	
ATA20	Locking Arms Refill, 20 pcs	
ATE12	Pouring/Ejector Bases Refill, 12 pcs	
ATK25	Articulation Keepers Refill, 25 pcs	
ATK300	Articulation Keepers Refill, 300 pcs	
ATS20	Articulation Spacers, 20 pcs	
ATT12	Acrylic frame Refill, 12 pcs	
ATT576	Acrylic frame Refill, 576 pcs	
ATN50	Transport System	
	1 × 50 pcs	Disposable Carrying Cases
	1 × 50 pcs	Plastic Retainers

PINDEX

Uncompromising model and dies precision assures parallel channels for quick, easy, and reliable pin placement. The laser drilling machine offers more accurate drill positioning for better die registration. The PINDEX Laser System reduces eye fatigue making easy to use with focus-free laser which outlasts conventional bulbs, and manual drilling allowing it to pre-set depth.



PINDEX Laser System



REF	Item
P113LCE	Laser System 230V Europe
	1 × Hand drill
	1 × Set of flexible socketformer
	2 × Tools for drill changing
	1 × Socket wrench
	1 × Pin introduction set
P113LUK	Laser System 230V UK
	1 × Hand drill
	1 × Set of flexible socketformer
	2 × Tools for drill changing
	1 × Socket wrench
	1 × Pin introduction set
PX113L	Laser System 115V
	1 × Hand drill
	1 × Set of flexible socketformer
	2 × Tools for drill changing
	1 × Socket wrench
	1 × Pin introduction set
PX125	Laser Retro-Fit Kit, 1 pcs



Drills



PX119

REF	Item
PX119	Carbide Drills, 3 pcs
PX120	Reamer, 1 pcs

Pins & Accessories



PX1125



PX146

REF	Item
PX112S	Medium Pins & Split Sleeves, 1000 pcs
PX112T	Medium Pins/Sleeve Brn, 1 pcs
PX114S	Long Pins, 1000 pcs
PX115	Short Pins & Plastic Sleeves, Bulk, 1000 pcs
PX116	Long Pin Kit, 100 pcs
PX117	Short Pins & Plastic Sleeves, 100 pcs
PX144	Self-Articulating Pins & Plastic Sleeves, 1000 pcs
PX145	Self-Articulating Pins & Plastic Sleeves, Bulk, 100 pcs
PX145S	Plastic Sleeves For Articulating Pins, Bulk, 1000 pcs
PX146	Flexible Rubber Base Molds (Full Arch, Right & Left Quadrants), 3 pcs
PX147	Flexible Rubber Base Molds (Full Arch), 1 pcs
PX148	Flexible Rubber Base Molds (Right/Left Quadrants), 2 pcs

Saw Blades




PX200

REF	Item
PX200	Handsaw Kit 1 × Handsaw 1 × 12 pcs Sawblades 0.18 mm 1 × 12 pcs Sawblades 0.25 mm
PX301	Saw Blades, .007" (.18mm), 100 pcs
PX302	Saw Blades, .010" (.25mm), 100 pcs
PX201	Saw Blades, .007" (.18mm), 12 pcs
PX202	Saw Blades, .010" (.25mm), 12 pcs
P38	PROXIMA SAW BLADES-12, 12 pcs

HYGENIC Waxes

Small details, big impact

The HYGENIC wax portfolio offers different wax sizes and shapes for various indications. Users can choose between bite wafers, u-shaped occlusal rim wax, square wax ropes, round wax strips, wax sticks which can be used for undercuts of individual trays, bite registration, bite plane appliances etc.



HYGENIC bite-wafers



H00825

This bite-mask contains a thin aluminium foil within two wax layers. U-shaped and designed for a fast precise retain of the situation.

REF	Item
H00825	HYGENIC Bite-Wafers, Double, 24 pcs
H00839	HYGENIC Bite-Wafers, Yellow, Soft, 50 pcs
H03365	HYGENIC Bite-Wafers, Green, Hard, 36 pcs

HYGENIC u-shaped occlusal rim wax



H04030

Recessed base aids adaption to base plate. Can be reshaped at room temperature to fit individual arch. If it's too rigid, just warm the wax slightly in warm water.

REF	Item
H00822	HYGENIC U-Shaped Occlusal Rim Wax, Standard Package, 100 pcs
H04030	HYGENIC U-Shaped Occlusal Rim Wax, Lab Package, 21 pcs

HYGENIC yellow bite wax



H00807

Economical alternative to beeswax with many of the same characteristics.

REF	Item
H00807	HYGENIC Yellow Bite Wax, 1 × 450 g
H00813	HYGENIC Yellow Bite Wax, 1 × 2.25 kg

HYGENIC boxing wax



H00815

Designed for boxing impressions, withstands extreme heat given off by dental stone during setting.

REF	Item
H00815	HYGENIC Boxing Wax, Regular, 1 × 450 g
H00816	HYGENIC Boxing Wax, Extra Fine, 1 × 450 g
H00837	HYGENIC Boxing Wax, Regular, 1 × 2.25 kg
H00838	HYGENIC Boxing Wax, Extra Fine, 1 × 2.25 kg

HYGENIC medium-soft wax



H00805

Easily softened and cooled to minimize shrinkage. Soft modeling.

REF	Item
H00806	HYGENIC Medium Soft, # 3 Pink Wax, 1 × 450 g
H00812	HYGENIC Medium Soft, # 3 Pink Wax, 1 × 2.25 kg
H00805	HYGENIC Red Wax, 1 × 450 g
H00811	HYGENIC Red Wax, 1 × 2.25 kg

HYGENIC orthodontic tray wax

Sticky, soft and narrow wax ropes.



H00827

REF	Item
H00827	HYGENIC Orthodontic Tray Wax, White, 48 pcs

HYGENIC wax ropes

Sticky, soft and narrow wax ropes.



H00819

REF	Item
H00819	HYGENIC Square Wax Ropes, Red, 1 pcs
H00817	HYGENIC Round Wax Strips, Red, 80 pcs
H00820	HYGENIC Square Wax Ropes, White, 1 pcs
H00818	HYGENIC Round Wax Strips, White, 80 pcs

HYGENIC wax sticks

Readily formed without heating, sticks are softer and tackier than „Round Wax Strips“ to minimize tissue displacement.



H00826

REF	Item
H00826	HYGENIC Wax Sticks, Scented, 48 pcs

HYGENIC extra tough wax

Durable, all-purpose wax that is a carefully prepared blend of exclusive waxes. Tough, pliable and sufficiently tacky to hold teeth well. Takes polish and boils out.



H00803

REF	Item
H00803	HYGENIC Extra Tough Wax, 1 x 450 g
H00809	HYGENIC Extra Tough Wax, 1 x 2.25 kg

Laboratory Accessories

COLTENEs laboratory accessories are made for laboratory needs. The HYGENIC Flexibole mixing bowls, as well as the nylon spatula are designed so that the mass to be mixed does not stick to them.



Laboratory Accessories



H00563

Plaster, stone, innvestments, and alginates will not stick to these self-cleaning Flexibole mixing bowls or the nylon spatula. More flexible than rubber bowls for more efficient mixing. The high-quality trimmer for silicones is used to remove any excess gingival mask material during the restoration process.



REF	Item
H00563	Mixing Bowl Small, 1 pcs
H00564	Mixing Bowl Medium, 1 pcs
H00565	Mixing Bowl Large, 1 pcs
H00566	Mixing Bowl Extra Large, 1 pcs
H00567	Mixing Bowl Jumbo, 1 pcs
H00676	3/16 x 3/32 x 50 Ft Amber Tube, 1 pcs
H00682	1/4 x 1/16 x 50 Ft Amber Tube, 1 pcs
H00684	1/4 x 1/16 x 50 Ft Black Tube, 1 pcs
H00685	1/4 x 3/32 Amber, 1 pcs
H00688	1/4 x 1/8 x 50 Ft Amber Tube, 1 pcs
H00697	5/16 X 3/32 X 50 Ft Amber Tube, 1 pcs
H00757	Saliva Eject Basketscreen A238, 1 pcs
H06012	Spatula, 3 pcs
H06541	100' Clean-Flex 1/2Hve Tubing, 1 pcs
8075	Trimmer for Silicones, 1 pcs

NOTES



LOCATIONS

CA	SciCan Ltd. 1440 Don Mills Road Toronto, Ontario, M3B 3P9 / Canada T +1 416 445 1600 F +1 416 445 2727	CH	Coltene/Whaledent AG Feldwiesenstrasse 20 9450 Altstätten / Switzerland T +41 71 757 53 00 F +41 71 757 53 01 info.ch@coltene.com
CN	Coltene/Whaledent China Trading CO. Ltd. B4046, Area B, 4th Floor, Building A No.9 Dongran North Street Beijing 100044 / P.R. China T +86 10 6212 3251 F +86 10 6212 3252 info.cn@coltene.com	DE	Coltene/Whaledent GmbH & Co. KG Raiffeisenstraße 30 89129 Langenau / Germany T +49 7345 805 0 F +49 7345 805 201 info.de@coltene.com
ES PT	Coltene Iberia S.L.U. c/Doctor Fleming, 3 - Bajo A 28220 Majadahonda, Madrid / Spain T +34 91 63 86 492 F +34 91 63 86 279 info.ib@coltene.com	FR	Coltene/Whaledent S.a.r.l. Synergie Park 3 avenue Louis Néel 59260 Lezennes / France T +33 320 64 64 00 F +33 320 64 64 03 info.fr@coltene.com
FR	Micro-Mega SA 12 rue du Tunnel 25000 Besancon / France T +33 381 54 42 42 F +33 381 54 42 39	IN	Coltene/Whaledent Private Ltd. Plot No. A2, MIDC, Pendhar Road, Taloja, Panvel Dist. Raigad. PIN-410208, Maharashtra / India T +91 22 2740 1841 info.in@coltene.com
IT	Coltene Italia S.r.l. c/o Polo Tecnologico della Brianza Via Alessandro Volta, 94 - Cond. Bianchi 20832 Desio (MB) / Italy T +39 0362 62 67 32 F +39 0362 63 25 80 info.it@coltene.com	JP	Coltene Japan LLC 2-2-28 2nd Floor, Nishikicho 190-0022 Tachikawa Tokyo / Japan T +81 42 595 6945 F +81 42 595 6946 info.jp@coltene.com
UK	Coltene/Whaledent Ltd. The President Suite Kendal House, Burgess Hill Victoria Way West Sussex, RH15 9NF / UK T +44 1444 235486 F +44 1444 870640 info.uk@coltene.com	US	Coltene/Whaledent Inc. 235 Ascot Parkway Cuyahoga Falls, Ohio 44223 / USA T +1 330 916 8800 F +1 330 916 7077 info.us@coltene.com



© COLTENE - www.coltene.com

